

Trichomonas Swab - TMA

Order Name: **Trich Swab**

Test Number: 6910121

Revision Date: 06/08/2018

TEST NAME	METHODOLOGY	LOINC CODE
Trichomonas Swab	Transcription-Mediated Amplification	46154-1

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	Swab	Swab	APTIMA Multitest Swab kit	Ambient / Refrigerated
Alternate 1	Swab	Swab	APTIMA Unisex Swab Collection kit	Ambient / Refrigerated
Instructions	Vaginal swab in APTIMA Multitest swab collection kit, or endocervical or male urethral swab in APTIMA Unisex Swab Specimen Collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A			

GENERAL INFORMATION

Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87661
Lab Section	Molecular Diagnostic