

## Mini-ANA Analyzer

Order Name: **Mini-ANA**  
Test Number: 6907869  
Revision Date: 10/01/2022

TEST NAME	METHODOLOGY	LOINC CODE
Mini-ANA Analyzer	Indirect Fluorescent Antibody	8061-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	4.5 mL (2 mL)	Serum	Clot Activator SST	Frozen
Alternate 1	4.5 mL (2 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Allow specimen to clot completely at room temperature. <b>Separate serum from cells ASAP or within 2 hours of collection and freeze for transport</b></p> <p>Stability After separation from cells: Ambient 8 hours, <b>Refrigerated 48 hours</b>, Frozen 3 months (avoid repeat freeze/thaw cycles).</p> <p><b>-- ! Updated Stability and Transport Temperature 9/3/2020 ! --</b></p>			

**GENERAL INFORMATION**

<b>Testing Schedule</b>	Mon - Fri
<b>Expected TAT</b>	3 - 7 Days
<b>Clinical Use</b>	To aid in the diagnosis of Connective Tissue Disease. This analyzer follows an algorithm or cascade of tests based on the results of the screening ANA test. The following tests will be ordered based on the ANA pattern(s) and titer: Anti-ds DNA, Anti-RNP, Anti-Sm, Anti-SSA Ro52, Anti-SSA Ro60, Anti-SSB/La, Anti-SCL-70, Jo-1.

<b>CPT Code(s)</b>	<h3>Initial Testing</h3>																						
	<table border="1"> <thead> <tr> <th>TEST NAME</th> <th>CPT CODES</th> </tr> </thead> <tbody> <tr> <td>ANA Screen</td> <td>86038</td> </tr> </tbody> </table>	TEST NAME	CPT CODES	ANA Screen	86038																		
TEST NAME	CPT CODES																						
ANA Screen	86038																						
	<h3>Possible Additional Testing</h3>																						
	<table border="1"> <thead> <tr> <th>TEST NAME</th> <th>CPT CODES</th> </tr> </thead> <tbody> <tr> <td>ANA Titer and Pattern</td> <td>86039</td> </tr> <tr> <td>Anti-Ribonuclear Antibody</td> <td>86235</td> </tr> <tr> <td>Anti-Smith Antibody</td> <td>86235</td> </tr> <tr> <td>Double Stranded DNA Antibody</td> <td>86225</td> </tr> <tr> <td>ENA Screen</td> <td>86235</td> </tr> <tr> <td>Jo-1 Antibody</td> <td>86235</td> </tr> <tr> <td>Scleroderma Antibody</td> <td>86235</td> </tr> <tr> <td>SSA/Ro52 Antibody (Sjogrens)</td> <td>86235</td> </tr> <tr> <td>SSA/Ro60 Antibody (Sjogrens)</td> <td>86235</td> </tr> <tr> <td>SSB/La Antibody (Sjogrens)</td> <td>86235</td> </tr> </tbody> </table>	TEST NAME	CPT CODES	ANA Titer and Pattern	86039	Anti-Ribonuclear Antibody	86235	Anti-Smith Antibody	86235	Double Stranded DNA Antibody	86225	ENA Screen	86235	Jo-1 Antibody	86235	Scleroderma Antibody	86235	SSA/Ro52 Antibody (Sjogrens)	86235	SSA/Ro60 Antibody (Sjogrens)	86235	SSB/La Antibody (Sjogrens)	86235
TEST NAME	CPT CODES																						
ANA Titer and Pattern	86039																						
Anti-Ribonuclear Antibody	86235																						
Anti-Smith Antibody	86235																						
Double Stranded DNA Antibody	86225																						
ENA Screen	86235																						
Jo-1 Antibody	86235																						
Scleroderma Antibody	86235																						
SSA/Ro52 Antibody (Sjogrens)	86235																						
SSA/Ro60 Antibody (Sjogrens)	86235																						
SSB/La Antibody (Sjogrens)	86235																						
<b>Lab Section</b>	Immunology - Serology																						