

## Candida and Trichomonas Vaginitis Swab

Order Name: **CV-TV Swab**  
 Test Number: 6987019  
 Revision Date: 04/06/2020

TEST NAME	METHODOLOGY	LOINC CODE
Candida Species - Swab	Transcription-Mediated Amplification	47000-5
Candida Glabrata - Swab	Transcription-Mediated Amplification	n/a
Trichomonas Swab (Vaginal)	Transcription-Mediated Amplification	70165-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
<b>Preferred</b>	<b>Swab</b>	<b>Swab</b>	<b>APTIMA Multitest Swab kit</b>	<b>Ambient / Refrigerated</b>
<b>Instructions</b>	Collect Swab in APTIMA Multitest Swab collection kit. Place swab in Swab Specimen Transport Tube, break shaft off at scoreline then recap tube. Stability: Room Temperature 30 Days, Refrigerated 30 Days, Frozen: N/A			

GENERAL INFORMATION	
Testing Schedule	Mon-Fri
Expected TAT	1-2 days after set-up
CPT Code(s)	87481, 87661
Lab Section	Molecular Diagnostic