Labcorp Oklahoma, Inc. Test Directory

Inheritest 14-Gene Panel

Order Name: Inherit 14
Test Number: 5194940
Revision Date: 03/21/2023

TEST NAME		METH	ODOLOGY	LOINC CODE	
Inheritest 14-Gene Panel		Polym	Polymerase Chain Reaction		
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	8.5 mL (3mL)	Whole Blood	ACD Solution A or B (Yellow Top)	Room Temperature	
Alternate 1	8.5 mL (3mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature	
Alternate 2	1	Saliva	Oragene Dx saliva kit	Room Temperature	
Alternate 3	1	Buccal swab	PurFlock buccal swab kit	Room Temperature	
Instructions	Specimen Type: Whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit Specimen Volume: 8.5 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit Mininum Volume: 3 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit Collection: Standard phlebotomy. Follow PurFlock buccal swab kit or Oragene Dx 500 saliva kit collection instructions. Do not eat, drink, smoke, or chew gum 30 min prior to collection. Specimen Storage: Maintain specimen at room temperature or refrigerate at 4C Do not freeze. Special Instructions: Males are not tested for x-linked disorders, including fragile X syndrome. Test orders must include an attestation that the				
Instructions	Specimen Volume: 8.5 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit Mininum Volume: 3 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit Collection: Standard phlebotomy. Follow PurFlock buccal swab kit or Oragene Dx 500 saliva kit collection instructions. Do not eat, drink, smoke, or chew gum 30 min prior to collection.				

GENERAL INFORMATION		
Expected TAT	14 - 21 days In some cases, additional time may be required for confirmatory or reflex tests.	
Performing Labcorp Test Code	481797	
Notes	Clinical Questionnaire for Inheritest® Carrier Screen and GeneSeq® PLUS	
Lab Section	Reference Lab	

provider has the patient's informed consent for genetic testing.

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025