## Labcorp Oklahoma, Inc. Test Directory

## Creatinine, Urine Random

Order Name: CREAT R U
Test Number: 3000750
Revision Date: 11/04/2025

| TEST NAME                |                       |   | METHODOLOGY                      | LOINC CODE            |  |
|--------------------------|-----------------------|---|----------------------------------|-----------------------|--|
| Creatinine, Urine Random |                       |   | Kinetic Alkaline Picrate (Jaffe) | 2161-8                |  |
| SPECIMEN REQUIREMENTS    |                       |   |                                  |                       |  |
| Specimen                 | Specimen Volume (min) | Specimen Type   | Specimen Container               | Transport Environment |  |
| Preferred                | 1 mL (0.5)            | Urine, Random   | Sterile Urine container          | Refrigerated          |  |
| Instructions             | ·                     | Random urine collection. No preservative. <b>Keep refrigerated.</b> Specimen stability: Ambient 3 days, Refrigerated 3 days, Frozen 2 weeks |                                  |                       |  |

| GENERAL INFORMATION |  |  |
|---------------------|--|--|
| Testing Schedule    | Daily  |  |
| Expected TAT        | 1-2 days   |  |
| Clinical Use        | Useful in monitoring renal function, to follow possible progression of renal disease, adjust dosages of medications in which renal excretion is pivotal. |  |
| CPT Code(s)         | 82570  |  |
| Lab Section         | Chemistry  |  |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025