

pH Serous Fluid

Order Name: **SRS PH**

Test Number: 3500350

Revision Date: 03/06/2020

| TEST NAME | METHODOLOGY | LOINC CODE |
|-----------------|-------------|------------|
| pH Serous Fluid | | 2748-2 |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|----------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.5) | Serous Fluid | Sterile Screwtop Container | On Ice |
| Instructions | Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Specimen must be on ice after collection. Deliver to lab immediately. | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Daily |
| Expected TAT | 1-2 days |
| Clinical Use | Useful in assessing acid-base balance. |
| CPT Code(s) | 83986 |
| Lab Section | Chemistry |