Alpha-fetoprotein (AFP) Serous Fluid

Order Name: SRS AFP
Test Number: 3620075
Revision Date: 06/12/2017

| TEST NAME | | | METHODOLOGY | LOINC CODE | |
|--------------------------------------|--|---------------|----------------------------|-----------------------|--|
| Alpha-fetoprotein (AFP) Serous Fluid | | | Chemiluminescence Assay | 11207-8 | |
| SPECIMEN REQUIREMENTS | | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 1 mL (0.5) | Serous Fluid | Sterile Screwtop Container | Refrigerated | |
| Instructions | Venous blood is often drawn simultaneously. Note fluid type on requisition and container. Stability: Ambient 8 hours. Refrigerated 48 hours. Freeze if testing cannot be started within 48 hours of collection. | | | | |

| GENERAL INFORMATION | | |
|---------------------|---|--|
| Testing Schedule | Mon, Wed, Fri | |
| Expected TAT | 1-3 days | |
| Clinical Use | Useful for the follow-up management of patients undergoing cancer therapy, especially for testicular and ovarian tumors and for primary hepatoma. | |
| CPT Code(s) | 82105 | |
| Lab Section | Chemistry | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025