

Amiodarone

Order Name: **AMIODARON**

Test Number: 3653300

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Amiodarone	Quant Liquid Chromatography-Tandem Mass Spectrometry	3330-8
N-Desethyl-Amiodarone	Quant Liquid Chromatography-Tandem Mass Spectrometry	6774-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Alternate 1	2 mL (0.5 mL)	Plasma	EDTA (Lavender Top)	Refrigerated
Instructions	<p>Specimen Type: Red-Top Tube, Lavender Top (EDTA) Tube OR Green-Top (Heparin) Tube.</p> <p>Specimen Storage: Submission/transport (<3 days): Room temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Specimen Collection: Transfer separated serum or plasma to a plastic transport tube before freezing.</p> <p>Specimen Stability: Ambient: 8 days, Refrigerated : 8 days, Frozen: 21 days</p>			

GENERAL INFORMATION	
Expected TAT	5 - 7 days
Clinical Use	Amiodarone is an antiarrhythmic drug. Therapeutic drug monitoring is useful to monitor compliance and avoid toxicity. Methodology: High Performance Liquid Chromatography (HPLC) or Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS).
Performing Labcorp Test Code	702891
Notes	Labcorp Test Code: 702891
CPT Code(s)	80151
Lab Section	Reference Lab