Labcorp Oklahoma, Inc. Test Directory

Tryptase

Order Name: **TRYPTASE**Test Number: 3658550
Revision Date: 01/06/2023

| TEST NAME | | | METHODOLOGY | LOINC CODE |
|----------------|---|---|--|-----------------------|
| Tryptase | | | ImmunoCAP | 21582-2 |
| SPECIMEN REQUI | REMENTS | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 1 mL (0.7 mL) | Serum | Clot Activator SST | Refrigerated |
| Alternate 1 | 1 mL (0.7 mL) | Serum | Clot Activator (Red Top, No-Gel) | Refrigerated |
| Instructions | Notes: 0.7 mL (Note: This volu Specimen Type: Gel-barrier, I Specimen Storage: Refrigera Specimen Collection: Separa Specimen Stability: Ambient: | lavender-top (EDTA) tube OR g te ate serum or plasma from cells | green-top (heparin) tube and transfer to a plastic transport tube. | |

| GENERAL INFORMATION | | |
|---------------------------------|---------------------------|--|
| Expected TAT | 2 - 4 days | |
| Performing Labcorp Test Code | 004280 | |
| Notes | Labcorp Test Code: 004280 | |
| CPT Code(s) | 83520 | |
| Lab Section | Reference Lab | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025