Neutrophil Membrane Antibody

Order Name: NEUTRO AB
Test Number: 5565160
Revision Date: 03/20/2023

| TEST NAME | 1 | | METHODOLOGY | LOINC CODE | |
|------------------------------|--|---------------|------------------------------------|-----------------------|--|
| Neutrophil Membrane Antibody | | | Flow cytometry | | |
| SPECIMEN REQUIREMENTS | | | | | |
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment | |
| Preferred | 2 mL (0.5) | Serum | Clot Activator (SST or Red No-Gel) | Frozen | |
| Instructions | Specimen Type: Gel-barrier tube or red-top tube Specimen Stability: Ambient: 5 days, Refrigerated : 5 days, Frozen: 1 year | | | | |

| GENERAL INFORMATION | | |
|---------------------------------|--|--|
| Expected TAT | 12-22 days | |
| Clinical Use | Evaluation for the presence of anti-neutrophil membrane antibody observed in some patients with neutropenia. | |
| Performing Labcorp Test Code | 832162 | |
| Notes | Labcorp Test Code: 832162 | |
| CPT Code(s) | 86021 | |
| Lab Section | Reference Lab | |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025