## **Acid Phosphate Total**

Order Name: ACID PHS T
Test Number: 2000280
Revision Date: 12/12/2022

TEST NAME		MET	THODOLOGY	LOINC CODE
Acid Phosphate Total		Enz	rymatic	12173-1
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.5 mL (0.8 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Allow specimen to clot completely. Transfer 1.5mL of serum into an aliquot tube and freeze.  CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.  Unacceptable Conditions: Plasma. Non-frozen specimens. Hemolyzed specimens.  Stability After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month			

GENERAL INFORMATION		
Expected TAT	5 - 8 days	
Clinical Use	This test is used in the evaluation of prostate cancer, Padget's disease, hyperthyroidism, and Gaucher's disease.	
Performing Labcorp Test Code	829642	
Notes	Labcorp Test Code: 829642	
CPT Code(s)	84060	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025