## **Tree Allergy Panel**

**Lab Section** 

Reference Lab

Order Name: AO TREE PN
Test Number: 2929545
Revision Date: 08/01/2023

TEST NAME			METHODOLOGY	LOINC CODE
Tree Allergy Panel			ImmunoCAP	
SPECIMEN REQUIRE	MENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5.6 mL (2.8 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			
GENERAL INFORMAT	TION 3-5 Days			
Notes	Box Elder (Maple) Tr Cedar/Red Tree IgE Common Silver Birch Cottonwood Tree All Elm Tree Allergen	Tree Allergen ergen	Hackberry IgE  Mesquite Tree Allergen  Mulberry Tree White Allergen  Oak Tree Allergen  Pecan, Hickory Tree Allergen	Sycamore, London Plane Tree Allergen  Walnut Tree Allergen  White Ash Tree Allergen  Willow Tree IgE
CPT Code(s)	86003 x14			

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025