

Penicillin Allergy Panel

Order Name: **PENICILIN**

Test Number: 5584380

Revision Date: 08/01/2023

TEST NAME	METHODOLOGY	LOINC CODE
Penicillin Allergy Panel	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1.6 mL (0.8 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION	
Expected TAT	3-5 Days
Performing Labcorp Test Code	805977, 805999, 602646, 067116
Notes	Testing Performed at Labcorp Phoenix as Individual Allergens.
CPT Code(s)	86003 x 4
Lab Section	Reference Lab