

Mini-ANA Analyzer

Order Name: **Mini-ANA**

Test Number: 6907869

Revision Date: 10/01/2022

TEST NAME	METHODOLOGY	LOINC CODE
Mini-ANA Analyzer	Indirect Fluorescent Antibody	8061-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	4.5 mL (2 mL)	Serum	Clot Activator SST	Frozen
Alternate 1	4.5 mL (2 mL)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	Allow specimen to clot completely at room temperature. <b>Separate serum from cells ASAP or within 2 hours of collection and freeze for transport</b> . Stability After separation from cells: Ambient 8 hours, <b>Refrigerated 48 hours</b> , Frozen 3 months (avoid repeat freeze/thaw cycles).  <b>-- ! Updated Stability and Transport Temperature 9/3/2020 ! --</b>			

GENERAL INFORMATION

Testing Schedule	Mon - Fri																										
Expected TAT	3 - 7 Days																										
Clinical Use	To aid in the diagnosis of Connective Tissue Disease. This analyzer follows an algorithm or cascade of tests based on the results of the screening ANA test. The following tests will be ordered based on the ANA pattern(s) and titer: Anti-ds DNA. Anti-RNP , Anti-Sm, Anti-SSA Ro52, Anti-SSA Ro60, Anti-SSB/La, Anti-SCL-70, Jo-1.																										
CPT Code(s)	<div><div>Initial Testing</div><table><tr><th>TEST NAME</th><th>CPT CODES</th></tr><tr><td>ANA Screen</td><td>86038</td></tr></table><div>Possible Additional Testing</div><table><tr><th>TEST NAME</th><th>CPT CODES</th></tr><tr><td>ANA Titer and Pattern</td><td>86039</td></tr><tr><td>Anti-Ribonuclear Antibody</td><td>86235</td></tr><tr><td>Anti-Smith Antibody</td><td>86235</td></tr><tr><td>Double Stranded DNA Antibody</td><td>86225</td></tr><tr><td>ENA Screen</td><td>86235</td></tr><tr><td>Jo-1 Antibody</td><td>86235</td></tr><tr><td>Scleroderma Antibody</td><td>86235</td></tr><tr><td>SSA/Ro52 Antibody (Sjogrens)</td><td>86235</td></tr><tr><td>SSA/Ro60 Antibody (Sjogrens)</td><td>86235</td></tr><tr><td>SSB/La Antibody (Sjogrens)</td><td>86235</td></tr></table></div>	TEST NAME	CPT CODES	ANA Screen	86038	TEST NAME	CPT CODES	ANA Titer and Pattern	86039	Anti-Ribonuclear Antibody	86235	Anti-Smith Antibody	86235	Double Stranded DNA Antibody	86225	ENA Screen	86235	Jo-1 Antibody	86235	Scleroderma Antibody	86235	SSA/Ro52 Antibody (Sjogrens)	86235	SSA/Ro60 Antibody (Sjogrens)	86235	SSB/La Antibody (Sjogrens)	86235
TEST NAME	CPT CODES																										
ANA Screen	86038																										
TEST NAME	CPT CODES																										
ANA Titer and Pattern	86039																										
Anti-Ribonuclear Antibody	86235																										
Anti-Smith Antibody	86235																										
Double Stranded DNA Antibody	86225																										
ENA Screen	86235																										
Jo-1 Antibody	86235																										
Scleroderma Antibody	86235																										
SSA/Ro52 Antibody (Sjogrens)	86235																										
SSA/Ro60 Antibody (Sjogrens)	86235																										
SSB/La Antibody (Sjogrens)	86235																										
Lab Section	Immunology - Serology																										