

Loxapine Level

Order Name: **LOXAPINE**

Test Number: 4001200

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Loxapine Level	Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Plasma	Lithium Heparin (Dark Green Top / No-Gel)	Refrigerated
Alternate 1	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Refrigerated
Instructions	<p><b>Specimen Type:</b> Red-top tube OR green-top (heparin) tube. <b>Gel-barrier tubes are NOT recommended.</b></p> <p><b>Specimen Storage:</b> SUBMISSION/TRANSPORT (&lt;3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p><b>Specimen Collection:</b> Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.</p> <p><b>Special Instructions:</b> Testing referred to Med Tox Laboratories Inc TC 65</p>			

GENERAL INFORMATION	
Expected TAT	5 - 10 days
Performing Labcorp Test Code	811240
Notes	Labcorp Test Code: 811240
CPT Code(s)	80299
Lab Section	Reference Lab