

## Eosinophil, Urine

Order Name: **EOS URINE**

Test Number: 1001600

Revision Date: 04/29/2019

TEST NAME	METHODOLOGY	LOINC CODE
Eosinophil, Urine	Wright Stain	12210-1

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Early morning specimens preferred. Keep specimen refrigerated, deliver to lab ASAP. Stability 24 Hours Refrigerated.			

### GENERAL INFORMATION

Testing Schedule	Daily
Expected TAT	1 Day
CPT Code(s)	85999
Lab Section	Hematology