Lobster Allergen

Order Name: LOBSTER
Test Number: 5610275
Revision Date: 02/11/2013

TEST NAME			METHODOLOGY	LOINC CODE	
Lobster Allergen			ImmunoCAP	15821-2	
SPECIMEN REQUIREMENTS					
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	0.4 mL (0.2 mL)	Serum	Clot Activator SST	Room Temperature	
Instructions	Specimen Type: Red-top tube or gel-barrier tube,				
	Separate Serum from Cells into a screwtop transport container.				
	Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)				

GENERAL INFORMATION		
Expected TAT	3-5 days	
CPT Code(s)	86003	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025