Labcorp Oklahoma, Inc. Test Directory

Penicillium Chrysogenym/notatum IgG

Order Name: **PEN CHRY G** Test Number: **5500433** Revision Date: **12/14/2022**

TEST NAME			METHODOLOGY	LOINC CODE
Penicillium Chrysogenym/notatum IgG			Enzyme immunoassay (EIA)	26957-
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.3 mL (0.1 mL)	Serum	Clot Activator SST	Room Temperature
Instructions	 Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy. Specimen Stability: Ambient: Not Available, Refrigerated : Not Available, Frozen: Not Available 			
GENERAL INFORMATION				

Expected TAT	2 - 4 days	
Performing Labcorp Test Code	605213	
Notes	Labcorp Test Code: 605213	
CPT Code(s)	86001	
Lab Section	Reference Lab	