

## Flow Cytometry on Body Fluid

Order Name: **FLOW FLUID**

Test Number: 8090085

Revision Date: 08/27/2024

TEST NAME	METHODOLOGY	LOINC CODE
Flow Cytometry on Body Fluid	Flow cytometry	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (3 mL)	Body Fluid	Sterile Screwtop Container	Refrigerated
Instructions	Aliquot CSF, Thoracentesis/Pleural fluid or Paracentesis/Ascites fluid in a sterile screwtop container. Keep Specimen Refrigerated. Indicate type of fluid on specimen tube and requisition. Please deliver to performing laboratory department (flow cytometry) ASAP. Specimen stability: 48hrs refrigerated <b>Rejection Criteria: WBC count - Minimum Requirement of 100 WBC/cmm. If less than 100 WBC/cmm, testing will be cancelled.</b>			

### GENERAL INFORMATION

Testing Schedule	Tue-Sat
Expected TAT	2 Days
CPT Code(s)	test dependant
Lab Section	Flow Cytometry