

Flow Cytometry on Tissue Biopsy

Order Name: **FLOW TISS**

Test Number: 8090065

Revision Date: 07/06/2022

TEST NAME	METHODOLOGY	LOINC CODE
Flow Cytometry on Tissue Biopsy	Flow cytometry	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	5x5mm	Tissue	Sterile Container with RPMI Solution	Refrigerated
Alternate 1	5x5mm	Tissue	Sterile Container with Saline	Refrigerated
Instructions	<p>Place a minimum of a 5x5mm tissue sample in a sterile screwtop container with RPMI. Keep specimen Refrigerated. Indicate type of tissue on specimen and requisition. Please deliver to performing laboratory department (flow cytometry) ASAP. DO NOT Place Tissue Into a Fixative Solution. Specimen stability: 48hrs refrigerated</p>			

GENERAL INFORMATION

Testing Schedule	Tue-Sat
Expected TAT	2 Days
CPT Code(s)	Test Dependent
Lab Section	Flow Cytometry