

Azithromycin Level by HPLC, Serum or Plasma

Order Name: **Azithromycin Lvl**

Test Number: 6904757

Revision Date: 01/02/2018

TEST NAME	METHODOLOGY	LOINC CODE
Azithromycin Level by HPLC, Serum or Plasma	Quantitative HPLC/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p>Specimen Type: Red-top tube Specimen Storage: FREEZE Specimen Collection: Serum should be separated from cells within two hours of venipuncture. Send serum in a plastic transport tube before freezing. Special Instructions: Testing referred to National Jewish Health NJHCO#PKAZM Was AZL NOTE: Drug dose and frequency, including time of last dose should be indicated on request form.</p>			

GENERAL INFORMATION	
Expected TAT	11-17 days
Performing Labcorp Test Code	857340
Notes	Labcorp Test Code: 857340
CPT Code(s)	80299
Lab Section	Reference Lab