

Lead, Whole Blood (Pediatric) Venous

Order Name: **Lead Peds Ven**
 Test Number: 6906338
 Revision Date: 12/09/2022

TEST NAME	METHODOLOGY	LOINC CODE
Lead, Whole Blood (Pediatric) Venous	Atomic Absorption	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Whole Blood	Royal Blue - EDTA (Trace Element Free)	Room Temperature
Alternate 1	1 mL (0.5 mL)	Whole Blood	EDTA K2 (Tan-Top tube)	Room Temperature
Instructions	<p>USE: Monitor environmental lead exposure in children younger than 16 years.</p> <p>Specimen Type: Royal blue-top (EDTA) tube OR Tan-top lead-free tube; submit original tube. Specimen Storage: Maintain specimen at room temperature. Specimen Collection: Mix blood thoroughly to avoid clotting. Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION

Expected TAT	1 - 4 days working day.
Performing Labcorp Test Code	717009
Notes	Labcorp Test Code: 717009
CPT Code(s)	83655
Lab Section	Reference Lab