

Serous Fluid, Routine Exam

Order Name: **SRS COUNT**

Test Number: 0800700

Revision Date: 02/26/2025

TEST NAME	METHODOLOGY	LOINC CODE
Color of Fluid	Visual	6824-7
Appearance of Fluid	Visual	9335-1
Total Nucleated Cells Body Fluid	Microscopy	n/a
WBC Count Body Fluid	Flow cytometry	26466-3
RBC Count Body Fluid	Flow cytometry	26455-6
Neutrophil Body Fluid Percent	Microscopy	26513-2
Lymphocyte Body Fluid Percent	Microscopy	11031-2
Eosinophil Body Fluid Percent	Microscopy	26452-3
Mononuclear Body Fluid Percent	Microscopy	26493-7
Id Pending Body Fluid	Visual	n/a
Specific Gravity Fluid	Optical	2964-5
Body Fluid Type	Prompt	47938-6

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3 mL (1 mL)	Body Fluid	EDTA (Lavender Top)	
Instructions	<p><b>Note fluid type on Tube and Request.</b> Good for all Body Fluids except for CSF, Semen and Synovial.</p> <p>Invert EDTA tube several times to mix well and deliver to laboratory immediately.</p> <p><b>Specimen stability</b> is 24hrs refrigerated preferred.</p> <p><b>Rejection Criteria:</b> Room Temperature specimens, Specimens received in syringes, urine cups, monovette tubes, red top tubes, gel-barrier tubes, or lithium heparin tubes; or specimens that are received frozen, clotted, or have exceeded stability.</p>			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
CPT Code(s)	89051
Lab Section	Urinalysis