

HCV FibroSure

Order Name: **HCV FibroSure**

Test Number: 5194382

Revision Date: 02/15/2023

TEST NAME	METHODOLOGY	LOINC CODE
HCV FibroSure	Assay Dependant	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3.5 mL (2 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated or Frozen
Instructions	<p><b>Specimen Type:</b> Gel-barrier tube is preferred. Red top tube is acceptable.</p> <p><b>Specimen Storage:</b> Specimen can be stored refrigerated at 2'C to 8'C for 72 hours and frozen at -70'C for 7 days. Frozen samples are stable for one freeze/thaw cycle.</p> <p><b>Specimen Collection:</b> Separate serum from cells within two hours of collection.</p> <p><b>Special Instructions:</b> The patient's age and gender must be included on the test request form.</p> <p><b>Specimen Stability:</b> Ambient: 72 hours, Refrigerated : 72 hours, Frozen: 7 days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 5 days
Performing Labcorp Test Code	550123
Notes	Labcorp Test Code: 550123
CPT Code(s)	81596
Lab Section	Reference Lab