

Drug Screen 11 w/Conf,Meconium

Order Name: **DRUG F S11**
Test Number: 5194626
Revision Date: 12/12/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|--------------------------------|------------------|------------|
| Drug Screen 11 w/Conf,Meconium | Immunoassay (IA) | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 5 g (1 g) | Meconium | Stool/Meconium Container | Refrigerated |
| Instructions | <p>Specimen Type: Meconium Collection Kit (PeopleSoft#96663) OR clean, leakproof polypropylene container.</p> <p>Specimen Storage: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.</p> <p>Specimen Collection: Specimens from different voidings may be pooled if necessary</p> <p>Special Instructions: Testing referred to MedTox Laboratories Inc MX#4077</p> | | | |

| GENERAL INFORMATION | |
|------------------------------|---------------|
| Expected TAT | 3 - 6 days |
| Performing Labcorp Test Code | 809558 |
| CPT Code(s) | 80307 |
| Lab Section | Reference Lab |